

MODEL	ITE 10-312	QUIX 10	RIXX RIC 10	MOOD RIC 312	MOOD Li-Ion RIC	SUN RIC 13	XS AHO 312	DUO AHO 13	DUO Li-Ion RIC	P AHO 13	HP AHO 675
Tech Level 16	TL 16 G5* ■▶*+Ⓜ	TL 16 G5 ■▶+Ⓜ	TL 16 G5 ▶+Ⓜ	TL 16 G5 ■▶*+Ⓜ	TL 16 G5 ■▶*+ⓂⓂ	TL 16 G5 ■▶*+Ⓜ		TL 16 G5 ■▶*+Ⓜ	TL 16 G5 ■▶*+ⓂⓂ	TL 16 G5 ■▶*+Ⓜ	TL 16 G4 ■●▶+Ⓜ
Tech Level 12	TL 12 G5* ■▶*+Ⓜ	TL 12 G5 ■▶+Ⓜ	TL 12 G5 ▶+Ⓜ	TL 12 G5 ■▶*+Ⓜ	TL 12 G5 ■▶*+ⓂⓂ	TL 12 G5 ■▶*+Ⓜ		TL 12 G5 ■▶*+Ⓜ	TL 12 G5 ■▶*+ⓂⓂ	TL 12 G5 ■▶*+Ⓜ	TL 12 G4 ■●▶+Ⓜ
Tech Level 8	TL 8 G5* ■▶*+Ⓜ	TL 8 G5 ■▶+Ⓜ	TL 8 G5 ▶+Ⓜ	TL 8 G5 ■▶*+Ⓜ	TL 8 G5 ■▶*+ⓂⓂ	TL 8 G5 ■▶*+Ⓜ	TL 8 G4 ●▶+Ⓜ	TL 8 G5 ■▶*+Ⓜ	TL 8 G5 ■▶*+ⓂⓂ	TL 8 G5 ■▶*+Ⓜ	TL 8 G4 ■●▶+Ⓜ
Tech Level 6	TL 6 G4 ●▶+Ⓜ	TL 6 G4 ●▶+Ⓜ	TL 6 G4 ▶+Ⓜ	TL 6 G4 ●▶+Ⓜ	TL 6 G5 ▶*+ⓂⓂ		TL 6 G4 ●▶+Ⓜ	TL 6 G4 ●▶+Ⓜ	TL 6 G5 ▶*+ⓂⓂ	TL 6 G4 ●▶+Ⓜ	TL 6 G4 ●▶+Ⓜ
Tech Level 4	TL 4 G4 ●▶+	TL 4 G4 ●▶+		TL 4 G4 ●▶+	TL 4 G5 ▶*+Ⓜ		TL 4 G4 ●▶+	TL 4 G4 ●▶+	TL 4 G5 ▶*+Ⓜ	TL 4 G4 ●▶+	TL 4 G4 ●▶+
Tech Level 3	TL 3 G4 ▶+	TL 3 G4 ▶+					TL 3 G4 ▶+	TL 3 G4 ▶+		TL 3 G4 ▶+	TL 3 G4 ▶+
Volta C		Volta C/PC ▶+					Volta C XS ▶+			Volta C P ▶+	Volta C HP ▶+

■

Kies de juiste CROS zender: (Bij ITE Wireless (WL) vereist)

CROS RIC G5 CROS Quix G5

CROS Li-Ion G5 CROS Li-Ion G4

●

Smart Connect, Smart Connect App, Smart Remote, Smart Transmitter compatible (Bij ITE Wireless (WL) vereist)

▶

Smart Key compatible

Ⓜ

Smart Direct App, Smart Mic, Smart Transmitter 2,4 compatible

+

Smart Remote App compatible

Ⓜ

Li-Ion Power vereist bij gebruik

- MOOD Li-Ion
- DUO Li-Ion
- CROS Li-Ion

🔊

Tinnitus

Beschikbaar in:

G4:

- Tinnitus noiser
- 5 statische ruissignalen
- 4 gemoduleerde ocean wave geluiden

G5:

- Tinnitus noiser
- 5 statische ruissignalen
- 4 gemoduleerde ocean wave geluiden
- Notch therapie

\* Alleen de IDA BT G5 is geschikt voor gebruik van Bluetooth accessoires en Smart Direct App.

MODÈLE	ITE 10-312	QUIX 10	RIXX RIC 10	MOOD RIC 312	MOOD Li-Ion RIC	SUN RIC 13	XS AHO 312	DUO AHO 13	DUO Li-Ion RIC	P AHO 13	HP AHO 675
Tech Level 16	TL 16 G5* ■▶*+Ⓜ	TL 16 G5 ■▶+Ⓜ	TL 16 G5 ▶+Ⓜ	TL 16 G5 ■▶*+Ⓜ	TL 16 G5 ■▶*+  Ⓜ	TL 16 G5 ■▶*+Ⓜ		TL 16 G5 ■▶*+Ⓜ	TL 16 G5 ■▶*+  Ⓜ	TL 16 G5 ■▶*+Ⓜ	TL 16 G4 ■●▶+Ⓜ
Tech Level 12	TL 12 G5* ■▶*+Ⓜ	TL 12 G5 ■▶+Ⓜ	TL 12 G5 ▶+Ⓜ	TL 12 G5 ■▶*+Ⓜ	TL 12 G5 ■▶*+  Ⓜ	TL 12 G5 ■▶*+Ⓜ		TL 12 G5 ■▶*+Ⓜ	TL 12 G5 ■▶*+  Ⓜ	TL 12 G5 ■▶*+Ⓜ	TL 12 G4 ■●▶+Ⓜ
Tech Level 8	TL 8 G5* ■▶*+Ⓜ	TL 8 G5 ■▶+Ⓜ	TL 8 G5 ▶+Ⓜ	TL 8 G5 ■▶*+Ⓜ	TL 8 G5 ■▶*+  Ⓜ	TL 8 G5 ■▶*+Ⓜ	TL 8 G4 ●▶+Ⓜ	TL 8 G5 ■▶*+Ⓜ	TL 8 G5 ■▶*+  Ⓜ	TL 8 G5 ■▶*+Ⓜ	TL 8 G4 ■●▶+Ⓜ
Tech Level 6	TL 6 G4 ●▶+Ⓜ	TL 6 G4 ●▶+Ⓜ	TL 6 G4 ▶+Ⓜ	TL 6 G4 ●▶+Ⓜ	TL 6 G5 ▶*+  Ⓜ		TL 6 G4 ●▶+Ⓜ	TL 6 G4 ●▶+Ⓜ	TL 6 G5 ▶*+  Ⓜ	TL 6 G4 ●▶+Ⓜ	TL 6 G4 ●▶+Ⓜ
Tech Level 4	TL 4 G4 ●▶+	TL 4 G4 ●▶+		TL 4 G4 ●▶+	TL 4 G5 ▶*+		TL 4 G4 ●▶+	TL 4 G4 ●▶+	TL 4 G5 ▶*+	TL 4 G4 ●▶+	TL 4 G4 ●▶+
Tech Level 3	TL 3 G4 ▶+	TL 3 G4 ▶+					TL 3 G4 ▶+	TL 3 G4 ▶+		TL 3 G4 ▶+	TL 3 G4 ▶+
Volta C		Volta C/PC ▶+					Volta C XS ▶+			Volta C P ▶+	Volta C HP ▶+

■  
Comment choisir le bon système CROS:  
(Au ITE Wireless (WL) nécessite)

CROS RIC G5 CROS Quix G5

CROS Li-Ion G5 CROS Li-Ion G4

●  
Smart Connect, Smart Connect App, Smart Remote, Smart Transmitter compatible (Au ITE Wireless (WL) nécessite)

▶  
Smart Key compatible

Ⓜ  
Smart Direct App, Smart Mic, Smart Transmitter 2,4 compatible

+

Smart Remote App compatible

||  
Li-Ion Power (puissance Li-Ion requise pour l'utilisation)

- MOOD Li-Ion
- DUO Li-Ion
- CROS Li-Ion

🔊  
Acouphène

Disponible en:

G4:

- Bruiteur anti-acouphènes
- 5 bruits statiques
- 4 bruits modulés de type vagues océaniques

G5:

- Bruiteur anti-acouphènes
- 5 bruits statiques
- 4 bruits modulés de type vagues océaniques
- Thérapie par Inhibition Latérale (TIL)

\* seul l'IDA BT G5 est compatible avec les accessoires Bluetooth et le Smart Direct App.