

MODEL	ITE 10-312	QUIX 10	RIXX RIC 10	MOOD RIC 312	MOOD Li-Ion RIC	SUN RIC 13	XS AHO 312	DUO AHO 13	DUO Li-Ion RIC	P AHO 13	HP AHO 675
											
Tech Level 16	TL 16 G5* ■ ▶ * + ◄	TL 16 G5 ■ ▶ + ◄	TL 16 G5 ▶ + ◄	TL 16 G5 ■ ▶ * + ◄	TL 16 G5 ■ ▶ * + II ◄	TL 16 G5 ■ ▶ * + ◄		TL 16 G5 ■ ▶ * + ◄	TL 16 G5 ■ ▶ * + II ◄	TL 16 G5 ■ ▶ * + ◄	TL 16 G4 ■ ● ▶ + ◄
Tech Level 12	TL 12 G5* ■ ▶ * + ◄	TL 12 G5 ■ ▶ + ◄	TL 12 G5 ▶ + ◄	TL 12 G5 ■ ▶ * + ◄	TL 12 G5 ■ ▶ * + II ◄	TL 12 G5 ■ ▶ * + ◄		TL 12 G5 ■ ▶ * + ◄	TL 12 G5 ■ ▶ * + II ◄	TL 12 G5 ■ ▶ * + ◄	TL 12 G4 ■ ● ▶ + ◄
Tech Level 8	TL 8 G5* ■ ▶ * + ◄	TL 8 G5 ■ ▶ + ◄	TL 8 G5 ▶ + ◄	TL 8 G5 ■ ▶ * + ◄	TL 8 G5 ■ ▶ * + II ◄	TL 8 G5 ■ ▶ * + ◄	TL 8 G4 ● ▶ + ◄	TL 8 G5 ■ ▶ * + ◄	TL 8 G5 ■ ▶ * + II ◄	TL 8 G5 ■ ▶ * + ◄	TL 8 G4 ■ ● ▶ + ◄
Tech Level 6	TL 6 G4 ● ▶ + ◄	TL 6 G4 ● ▶ +	TL 6 G4 ▶ + ◄	TL 6 G4 ● ▶ + ◄	TL 6 G5 ▶ * + II ◄		TL 6 G4 ● ▶ + ◄	TL 6 G4 ● ▶ + ◄	TL 6 G5 ▶ * + II ◄	TL 6 G4 ● ▶ + ◄	TL 6 G4 ● ▶ + ◄
Tech Level 4	TL 4 G4 ● ▶ +	TL 4 G4 ● ▶ +		TL 4 G4 ● ▶ +	TL 4 G5 ▶ * + II		TL 4 G4 ● ▶ +	TL 4 G4 ● ▶ +	TL 4 G5 ▶ * + II	TL 4 G4 ● ▶ +	TL 4 G4 ● ▶ +
Tech Level 3	TL 3 G4 ▶ +	TL 3 G4 ▶ +					TL 3 G4 ▶ +	TL 3 G4 ▶ +		TL 3 G4 ▶ +	TL 3 G4 ▶ +
Volta C		Volta C/PC ▶ +					Volta C XS ▶ +			Volta C P ▶ +	Volta C HP ▶ +

■ Kies de juiste CROS zender:
(Bij ITE Wireless (WL) vereist)




CROS RIC G5 CROS Li-ion G5

● Smart Connect,
Smart Connect App, Smart Remote, Smart Transmitter compatible (Bij ITE Wireless (WL) vereist)






CROS Quix G5

▶ Smart Key compatible



luetooth

Smart Direct App, Smart Mic, Smart Transmitter 2,4 compatible




+

Smart Remote App compatible



||

Li-Ion Power vereist bij gebruik

- MOOD Li-Ion
- DUO Li-Ion
- CROS Li-Ion



🔊

Tinnitus

Beschikbaar in:

G4:

- Tinnitus noiser
- 5 statische ruissignalen
- 4 gemoduleerde ocean wave geluiden

G5:

- Tinnitus noiser
- 5 statische ruissignalen
- 4 gemoduleerde ocean wave geluiden
- Notch therapie

MODÉLE	ITE 10-312	QUIX 10	RIXX RIC 10	MOOD RIC 312	MOOD Li-Ion RIC	SUN RIC 13	XS AHO 312	DUO AHO 13	DUO Li-Ion RIC	P AHO 13	HP AHO 675
											
Tech Level 16	TL 16 G5* ■ ▶ * + ◄	TL 16 G5 ■ ▶ + ◄	TL 16 G5 ▶ + ◄	TL 16 G5 ■ ▶ * + ◄	TL 16 G5 ■ ▶ * + II ◄	TL 16 G5 ■ ▶ * + ◄		TL 16 G5 ■ ▶ * + ◄	TL 16 G5 ■ ▶ * + II ◄	TL 16 G5 ■ ▶ * + ◄	TL 16 G4 ■ ▷ + ◄
Tech Level 12	TL 12 G5* ■ ▶ * + ◄	TL 12 G5 ■ ▶ + ◄	TL 12 G5 ▶ + ◄	TL 12 G5 ■ ▶ * + ◄	TL 12 G5 ■ ▶ * + II ◄	TL 12 G5 ■ ▶ * + ◄		TL 12 G5 ■ ▶ * + ◄	TL 12 G5 ■ ▶ * + II ◄	TL 12 G5 ■ ▶ * + ◄	TL 12 G4 ■ ▷ + ◄
Tech Level 8	TL 8 G5* ■ ▶ * + ◄	TL 8 G5 ■ ▶ + ◄	TL 8 G5 ▶ + ◄	TL 8 G5 ■ ▶ * + ◄	TL 8 G5 ■ ▶ * + II ◄	TL 8 G5 ■ ▶ * + ◄	TL 8 G4 ● ▶ + ◄	TL 8 G5 ■ ▶ * + ◄	TL 8 G5 ■ ▶ * + II ◄	TL 8 G5 ■ ▶ * + ◄	TL 8 G4 ■ ▷ + ◄
Tech Level 6	TL 6 G4 ● ▷ + ◄	TL 6 G4 ● ▷ + ◄	TL 6 G4 ▶ + ◄	TL 6 G4 ● ▷ + ◄	TL 6 G5 ▶ * + II ◄		TL 6 G4 ● ▷ + ◄	TL 6 G4 ● ▷ + ◄	TL 6 G5 ▶ * + II ◄	TL 6 G4 ● ▷ + ◄	TL 6 G4 ● ▷ + ◄
Tech Level 4	TL 4 G4 ● ▷ +	TL 4 G4 ● ▷ +		TL 4 G4 ● ▷ +	TL 4 G5 ▶ * + II		TL 4 G4 ● ▷ +	TL 4 G4 ● ▷ +	TL 4 G5 ▶ * + II	TL 4 G4 ● ▷ +	TL 4 G4 ● ▷ +
Tech Level 3	TL 3 G4 ▶ +	TL 3 G4 ▶ +					TL 3 G4 ▶ +	TL 3 G4 ▶ +		TL 3 G4 ▶ +	TL 3 G4 ▶ +
Volta C		Volta C/PC ▶ +					Volta C XS ▶ +			Volta C P ▶ +	Volta C HP ▶ +

■ Comment choisir le bon système CROS:
(Au ITE Wireless (WL) nécessite)




CROS RIC G5 CROS Li-ion G5

● Smart Connect,
Smart Connect App, Smart
Remote, Smart Transmitter
compatible (Au ITE Wireless
(WL) nécessite)






CROS Quix G5

▶ Smart Key compatible



luetooth

Smart Direct App, Smart
Mic, Smart Transmitter 2,4
compatible




+

Smart Remote App compatible



||

Li-Ion Power
(puissance Li-Ion requise pour l'utilisation)

- MOOD Li-Ion
- DUO Li-Ion
- CROS Li-Ion



🔊

Acouphène

Disponible en:

G4:

- Bruiteur anti-acouphènes
- 5 bruits statiques
- 4 bruits modulés de type vagues océaniques

G5:

- Bruiteur anti-acouphènes
- 5 bruits statiques
- 4 bruits modulés de type vagues océaniques
- Thérapie par Inhibition Latérale (TIL)